



# Second Skin

CUSTOMISED COMPRESSION GARMENTS & SPLINTS


## Washing tips

### SOILED GARMENTS & SPLINTS (BLOOD STAINS\*)



- For fresh blood stains simply **work in a small amount of heavy-duty liquid laundry detergent to break apart the stain**. Work the detergent into the stain with your fingers or a soft bristled brush.

(For dried or old blood stains, scrape off any crusted matter with a spoon)



- **Allow the detergent to work into the stain for 15 minutes**, then wash as recommended on the fabric care label.

- Check the stained area after washing, **do not dry the area with a heated appliance** (e.g. hairdryer, as heat will lock in the stain). If traces of the stain remain move to the next step below.



- **Mix a solution of cool water and oxygen-based stain remover** (e.g. Vanish or Napisan, following the directions to use). Completely submerge the blood-stained item, soak for at least two hours.



- **Wash as usual in warm or cold water with detergent**. If the stain remains, treat again by soaking in a fresh solution, then rewash.



- Once the garment or splint is clean, **dry flat in an open ventilated area**

\* Please take the necessary precautions, such as wearing gloves, when coming into contact with bodily fluids